

FULL PROTECTION

APRON FEATURES

- Adjustable Buckle Set**
Adjustable buckle closures ensure a secure fit.
- Back Relief Belt**
Elastic lumbar support belt with buckle closure system.
- Seam-to-seam Protection**
Seam-to-seam, standardized protection in Front (no change in protection at overlap)

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro and Buckle Closures

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Wings	Extra Pocket
Embroidery	Logo	Sleeve

CORE MATERIAL

THICKNESS OPTIONS

 Standard Lead	0.25, 0.50mm, 1.00mm
 Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
 Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
 Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
 Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

BALANCE™ WRAPAROUND BACK RELIEF APRON

MODEL: [AP-WABR]

Offering full protection, this back relieving design transfers the weight from your shoulders to your hips with the aid of a 6" lumbar support belt. The Wrap Around Back Relief apron features velcro and buckle closures to ensure a proper fit. The Essential apron configuration **features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection. The Essential Wraparound is available in standard sizing only.** For a made-to-measure Wraparound apron, see Balance Wraparound.



Example of a 0.50mm LE
Balance Vest Skirt Apron
Lead Configuration

