

FRONTAL PROTECTION

APRON FEATURES

+ Anti-slip fabric (interior)
Reduces slippage and improves fit.

+ Built-in Thyroid Collar
A built-in, velcro thyroid collar provides added protection for your thyroid

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

No Closure

Available in adult and pediatric sizes

ADDITIONAL CUSTOMIZATION OPTIONS

Embroidery Logo

CORE MATERIAL

THICKNESS OPTIONS

Pb	Standard Lead	0.25, 0.50mm, 1.00mm
L	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
LF	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
PR	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
B	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

DENTAL APRON

MODEL: [AP-DA]

Features a non-slip microfiber material that ensures the apron lays in place over the patient. There are no closures on the dental apron. A built-in thyroid collar is included.

Adult sizing: 28"W x 31"L; Pediatric sizing: 24"W x 23"L

