

FULL PROTECTION

APRON FEATURES

- +** **Adjustable Buckle Set**
Adjustable buckle closures ensure a secure fit.
- +** **Back Relief Belt**
Elastic lumbar support belt with buckle closure system.
- +** **Extended Front Panels**
Front panels are extended to 3/4 overlap coverage where LE protection is doubled

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro and Buckle Closures



ESSENTIAL™ WRAPAROUND BACK RELIEF APRON

MODEL: [AP-WAFBMA]

Offering full protection, this back relieving design transfers the weight from your shoulders to your hips with the aid of a 6" lumbar support belt. The Wrap Around Back Relief apron features velcro and buckle closures to ensure a proper fit. The Essential apron configuration **features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection. The Essential Wraparound is available in standard sizing only.** For a made-to-measure Wraparound apron, see Balance Wraparound.

ADDITIONAL CUSTOMIZATION OPTIONS

- | | | |
|----------------|-------|--------------|
| Thyroid Collar | Wings | Extra Pocket |
| Embroidery | Logo | Sleeve |

CORE MATERIAL

THICKNESS OPTIONS

Pb	Standard Lead	0.25, 0.50mm, 1.00mm
U	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
P	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
PR	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
B	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.



Example of a 0.50mm LE Essential Vest Skirt Apron Lead Configuration

