

FRONT PROTECTION

APRON FEATURES



Backstrap

Velcro backstrap securely fastens shoulder area.



Elastic Back

Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.

Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar

Drop Shoulder

Wings

Extra Pocket

Embroidery

Back Mesh

Logo

Sleeve



FLEXBACK APRON

MODEL: [AP-BR]

Representing the gold standard of quality, the Flexback apron design provides enhanced back relief, comfort, and excellent weight distribution. Criss-cross back straps are specially designed with elastic to “flex” the apron’s weight burden off your shoulders to your hips, thereby providing support to your lower back. A velcro backstrap (shown below) is included and secures the apron around your shoulder area to prevent sliding.

CORE MATERIAL

THICKNESS OPTIONS

	Standard Lead	0.25, 0.50mm, 1.00mm
	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

