

FULL PROTECTION

## APRON FEATURES

- +** **Elastic Back**  
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.
- +** **Single Front Panel**  
Eliminates excess overlapping panels achieving a 10% weight reduction (compared to VS model).
- +** **Tri-tab Closure Design**  
Ensures a proper fit.
- +** **Non-slip Backing (skirt)**  
Reduces slippage and improves fit.

- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Velcro Closure**

### ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Extra Pocket
- Embroidery
- Logo
- Sleeve

### CORE MATERIAL

### THICKNESS OPTIONS

<b>Pb</b>	Standard Lead	0.25, 0.50mm, 1.00mm
<b>L</b>	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
<b>LF</b>	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
<b>PR</b>	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
<b>B</b>	Bilayer	0.25, 0.50mm, 1.00mm



# TRI-TAB VEST SKIRT APRON

MODEL: [ AP-TVS ]

With an easy-on, easy-off multi tab (velcro) closure design, the Tri-tab is a popular style for vest and skirt combinations. Velcro closures are located on the left shoulder and around the hips. An elastic back allows for a flexible fit while providing premium comfort and back support. The Tri-tab Vest Skirt reduces excess overlapping lead and is **up to 10% lighter than traditional vest models.**



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.