

FRONT PROTECTION

APRON FEATURES



Wide Body Back Panels
The comfortable panels secure your apron and provides lumbar support.

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket


Velcro Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Drop Shoulder	Wings
Extra Pocket	Embroidery	Backstrap
Logo	Sleeve	Back Mesh

CORE MATERIAL

THICKNESS OPTIONS

 Standard Lead	0.25, 0.50mm, 1.00mm
 Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
 Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
 Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
 Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

VELCRO ADJUSTABLE APRON

MODEL: [AP-VA]

Do you need comfort and adjustability at an affordable price? The Velcro Adjustable apron has crisscrossing back straps that provide back support while securing the apron around your waist. Wide back straps prevent the apron from slipping around your shoulder and back area. Pediatric sizing available.

