

FULL PROTECTION

## APRON FEATURES

- +** **Elastic Back**  
Provides enhanced adjustability and lumbar support. Shifts weight from shoulders & hips.
- +** **Single Front Panel**  
Eliminates excess overlapping panels achieving a 15% weight reduction (compared to other wraparound and VS models).
- +** **Tri-tab Closure**  
Ensures a proper fit.

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Velcro Closure**

### ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Wings	Extra Pocket
Embroidery	Logo	Sleeve



# WRAPAROUND FLEXBACK APRON

MODEL: [ AP-WAFB ]

Want the ultimate in comfort and protection? The Wraparound Flexback apron combines the best qualities of our flexback apron with the full protection and back support of our wraparound back relief apron. The elastic stretch-back shifts weight from your shoulders to your hips while the tri-tab velcro closure system provides convenient, easy-on easy-off wearability. The Wrap Around Flexback eliminates excess overlapping lead and is up to **15% lighter than traditional wrap styles.**

### CORE MATERIAL

### THICKNESS OPTIONS

<b>Pb</b>	Standard Lead	0.25, 0.50mm, 1.00mm
<b>L</b>	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
<b>PR</b>	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
<b>PR</b>	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
<b>B</b>	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.