

FRONT PROTECTION

APRON FEATURES



Backstrap

Velcro backstrap securely fastens shoulder area.



Back Relief Belt

Wide, elastic belt provides lumbar support and maximum comfort.

Shoulder Padding

Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Belt Buckle Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar

Drop Shoulder

Wings

Extra Pocket

Embroidery

Back Mesh

Logo

Sleeve

CORE MATERIAL

THICKNESS OPTIONS

	Standard Lead	0.25, 0.50mm, 1.00mm
	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

BACK RELIEF APRON

MODEL: [AP-BR]

Save yourself from back pain with Protech's Back Relief apron. This apron features a unique padded belt system constructed with a bonded elasticized panel that will not lose its shape. When engaged, the special back relief belt provides lumbar support while transferring the weight of the apron from your shoulders to your hips for increased comfort and back support throughout the life of the apron. A velcro backstrap is also included to secure the apron around your shoulder area and prevent sliding.

