

FULL PROTECTION

APRON FEATURES



- Adjustable Buckle Set**
Adjustable buckle closures ensure a secure fit.
- Non-slip Backing (skirt interior)**
Reduces slippage and improves fit
- Seam-to-seam Protection**
Seam-to-seam, standardized protection in Front (no change in protection at overlap)

- Shoulder Padding**
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket**
- Velcro and Buckle Closures**

ADDITIONAL CUSTOMIZATION OPTIONS

- Thyroid Collar
- Wings
- Extra Pocket
- Embroidery
- Logo
- Sleeve

CORE MATERIAL

THICKNESS OPTIONS

Pb	Standard Lead	0.25, 0.50mm, 1.00mm
L	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
P	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
PR	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
B	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

BALANCE™ VEST SKIRT APRON

MODEL: [AP-VS]

Be prepared for prolonged procedures with the Vest Skirt apron. The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue. The Balance configuration features **seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration). The skirt comes standard with a premium non-slip fabric in the back.



Example of a 0.35mm LE Balance Vest Skirt Apron Lead Configuration

