

FULL PROTECTION

APRON FEATURES



Adjustable Buckle Set
Adjustable buckle closures ensure a secure fit.



Flexback Belt
Wraps internally around your back and securely fastens in the front.



Non-slip Backing (skirt interior)
Reduces slippage and improves fit.



Seam-to-seam Protection
Seam-to-seam, standardized protection in Front (no change in protection at overlap)

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro and Buckle Closures

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar
Embroidery

Wings
Logo

Extra Pocket
Sleeve

CORE MATERIAL

THICKNESS OPTIONS

	Standard Lead	0.25, 0.50mm, 1.00mm
	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
	Bilayer	0.25, 0.50mm, 1.00mm




BALANCE™ VEST SKIRT FLEXBACK APRON

MODEL: [AP-VSFB]

The Vest Skirt Flexback shifts 70% of the weight directly to the hips with added support of an adjustable, elastic back relief belt. The Balance configuration **features seam-to-seam, single LE protection with standard overlapping panels** (see below for lead configuration). The skirt comes standard with an anti-slip fabric in the back (interior) to keep your skirt in place.

Example of a 0.50mm LE Balance Vest Skirt Apron Lead Configuration



0.50mm LE

0.50mm LE

0.25mm LE



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.