

## FRONT PROTECTION

# APRON FEATURES



**Wide Body Back Panels**  
The comfortable panels secure your apron and provides lumbar support.

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Durable Buckle Closure**

# BUCKLE APRON

MODEL: [ AP-BA ]

The Buckle Apron comes standard with an easy-release buckle allowing for quick removal of the apron. The belt is adjustable to accommodate different sizes. Criss-cross back straps provide a secure, comfortable fit.

## ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Drop Shoulder	Elastic Back
Extra Pocket	Embroidery	Backstrap
Logo	Sleeve	

## CORE MATERIAL

## THICKNESS OPTIONS

 Standard Lead	0.25, 0.50mm, 1.00mm
 Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
 Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
 Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
 Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

