

FULL PROTECTION

APRON FEATURES

- +** Adjustable Buckle Set
Adjustable buckle closures ensure a secure fit.
- +** Non-slip Backing (skirt interior)
Reduces slippage and improves fit
- +** Extended Front Panels
Front panels are extended to 3/4 overlap coverage where LE protection is doubled

- Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.
- Apron Pocket
- Velcro and Buckle Closures

ADDITIONAL CUSTOMIZATION OPTIONS

- | | | |
|----------------|-------|--------------|
| Thyroid Collar | Wings | Extra Pocket |
| Embroidery | Logo | Sleeve |

CORE MATERIAL

THICKNESS OPTIONS

Pb	Standard Lead	0.25, 0.50mm, 1.00mm
L	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
LF	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
PR	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
B	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

ESSENTIAL™ VEST SKIRT APRON

MODEL: [AP-VSQ]

Be prepared for prolonged procedures with the Vest Skirt apron. The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue. The Essential apron configuration features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection. The skirt comes standard with a premium non-slip fabric in the back. Essential aprons are only available in standard sizes. For a made-to-measure Vest Skirt, see the Balance Vest Skirt Apron.



Example of a 0.35mm LE Essential Vest Skirt Apron Lead Configuration

