

FULL PROTECTION

## APRON FEATURES

- +** **Adjustable Buckle Set**  
Adjustable buckle closures ensure a secure fit.
  - +** **Flexback Belt**  
Wraps internally around your back and securely fastens in the front.
  - +** **Non-slip Backing (skirt interior)**  
Reduces slippage and improves fit.
  - +** **Extended Front Panels**  
Front panels are extended to 3/4 overlap coverage where LE protection is doubled
- Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.
  - Apron Pocket**
  - Velcro and Buckle Closures**



# ESSENTIAL™ VEST SKIRT FLEXBACK APRON

MODEL: [ AP-VSFBQ ]

The Essential Vest Skirt Flexback shifts 70% of the weight directly to the hips with added support of an adjustable, elastic back relief belt. The Essential apron configuration features **extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection.** The skirt comes standard with an anti-slip fabric in the back (interior) to keep your skirt in place. **Essential aprons are only available in standard sizes.** For a made-to-measure Vest Skirt Flexback, see the Balance Vest Skirt Flexback Apron.

### ADDITIONAL CUSTOMIZATION OPTIONS

- |                |       |              |
|----------------|-------|--------------|
| Thyroid Collar | Wings | Extra Pocket |
| Embroidery     | Logo  | Sleeve       |

### CORE MATERIAL

### THICKNESS OPTIONS

<b>Pb</b>	Standard Lead	0.25, 0.50mm, 1.00mm
<b>U</b>	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
<b>P</b>	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
<b>PR</b>	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
<b>B</b>	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.



### Example of a 0.50mm LE Essential Vest Skirt Apron Lead Configuration

