

FULL PROTECTION

APRON FEATURES

+ Complete Coverage
The front & back panels ensure you are fully protected from beginning to end of your pregnancy.

+ Half Apron Included
A Half apron is included for additional protection as you progress through your pregnancy.

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Tie & Buckle Closures

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar Wings
Extra Pocket Embroidery
Logo

CORE MATERIAL

THICKNESS OPTIONS

	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
	Bilayer	0.25, 0.50mm, 1.00mm

MATERNITY APRON

MODEL: [AP-MA]

The maternity apron features easy ties and a quick release buckle to secure the back and front panels; allowing complete adjustability during pregnancy. The maternity apron also includes an additional .50mm equivalent half apron for a total of 1.00mm Pb protection over the abdominal area and lower extremities.



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.