

FRONT PROTECTION

APRON FEATURES

+ **Hanging Straps**
Hanging straps on the shoulders allow for easy storage on any Protech apron rack.

+ **Surgical Drop**
Quick release velcro closure tabs help maintain sterile field.

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro Closure

MODEL: [AP-SD]

Maintain your sterile field with the quick release Surgical Drop Apron. Easy-release velcro closures allow the surgical drop to fall even when underneath your scrubs. With padded shoulders and a criss-cross back strap design the surgical drop apron provides maximum comfort for all users. Hanging straps are included to allow for easy storage on any Protech apron rack.

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar Extra Pocket Embroidery
Logo

CORE MATERIAL

THICKNESS OPTIONS

Pb	Standard Lead	0.25, 0.50mm, 1.00mm
L	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
LF	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
PR	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
B	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

