

## FRONT PROTECTION

# APRON FEATURES

**Shoulder Padding**  
Premium foam padding in the shoulders provides extra comfort.

**Apron Pocket**

**Tie Closure**

## ADDITIONAL CUSTOMIZATION OPTIONS


Thyroid Collar  
Logo

Embroidery  
Sleeve

Wings  
Backstrap

## CORE MATERIAL

## THICKNESS OPTIONS

 Standard Lead	0.25, 0.50mm, 1.00mm
 Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
 Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
 Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
 Bilayer	0.25, 0.50mm, 1.00mm

The Tie Apron is an easy-to-use, affordable option for x-ray protection. Utilizing a tie closure, the Tie Apron's criss crossing straps ensure a secure and comfortable fit. Pediatric sizing available, see page 52.

# TIE APRON

MODEL: [ AP-TA ]



For sizing and color options, please see our sizing and fabric color documents or view them online at [protechmed.com](http://protechmed.com).



\*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.