

FRONT PROTECTION

APRON FEATURES

+ Multi-Panel Design
Allows for easy maneuvering and a comfortable fit.

+ Backstrap
Velcro backstrap securely fastens shoulder area.

Shoulder Padding
Premium foam padding in the shoulders provides extra comfort.

Apron Pocket

Velcro Closure

ADDITIONAL CUSTOMIZATION OPTIONS

Thyroid Collar	Drop Shoulder	Wings
Extra Pocket	Embroidery	Back Mesh
Logo	Sleeve	

If you sit while performing procedures requiring x-ray, you may be unnecessarily exposing yourself to radiation. The Urology Apron features a three panel system at the bottom of the apron which enables the user to sit comfortably while maintaining maximum coverage from x-ray.

CORE MATERIAL

THICKNESS OPTIONS

Pb	Standard Lead	0.25, 0.50mm, 1.00mm
L	Ultralite Lead	0.25, 0.35, 0.50mm, 1.00mm
PR	Prolite Lead-Free	0.25, 0.35, 0.50mm, 1.00mm
PR	Premier Lead Free	0.25, 0.35, 0.50mm, 1.00mm
B	Bilayer	0.25, 0.50mm, 1.00mm



For sizing and color options, please see our sizing and fabric color documents or view them online at protechmed.com.



*CE applies to Lead, Ultralite Lead and Prolite Lead Free IEC in 0.25, 0.35, and 0.50mm. Prolite IEC is limited to 110kVp diagnostic range.

UROLOGY APRON

MODEL: [AP-UA]

