



**PROTECH**  
Medical

**PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# 42 FITOVER

MODEL: [ RE-42BLK ]



## DESCRIPTION:

The 42 fitover is made to wear over prescription eyewear that are wider or more rectangular in shape. It features a smaller vertical dimension than other fitover styles as well as 0.5mm Pb lead sheeting for side shields, thus decreasing the overall weight by approximately 20%.

## DESIGN:

**Frame style:** Fitover

**Front lens:** 0.75mm Pb Lead Glass

**Frame Material:** Plastic

**Weight:** 92g

**Frame size:** 63-14-135

**Side Shields:** 0.50mm (standard) Pb Strips

## PRESCRIPTION:

**Options:** N/A

**Range:** N/A

**Minimum Progressive PD:** N/A

**Base Curve:** 4

## FEATURES:

Large side shields included

## COLORS:

Black (BLK)



*Side Shields (0.50mm LE) Included*