



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

53 WRAP

MODEL: [RE-53 / 53A]



DESCRIPTION:

The 53 wrap is constructed of durable nylon that exhibits superior flexibility while minimizing potential breakage. Its full saddle bridge provides maximum comfort and evenly distributes weight over the nose. Also available with adjustable nose pads (Model 53A).

DESIGN:

Frame style: Classic
Front lens: 0.75mm Pb Glass
Frame Material: Nylon Plastic

Weight: 84g
Frame size: 58-16-135
Side Shields: 0.50mm (standard) or 0.75mm Pb Lead Glass

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -10.00

Minimum Progressive PD: 62
Base Curve: 4

FEATURES:

Adjustable Nose Pads, Spring Hinge Temples

COLORS:

Black (BLK), Tortoise (TORT)

