



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# 553S METALITE

MODEL: [ RE-553S ]

## DESCRIPTION:

The model 553 is a metal frame featuring custom fitted 0.50mm lead strip side shields that create a seal to the front lenses for complete lateral coverage. Its silver color and retro style offer a classic look while non-slip, silicone nose pads provide easy adjustability and allow the frame to fit a wide variety of face types.

## DESIGN:

Frame style: Metals

Front lens: 0.75mm Pb Glass

Frame Material: Optical Steel

Weight: 62g

Frame size: 50-18-140

Side Shields: 0.50mm (included)

## PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive

Range: +8.00 to -9.00

Minimum Progressive PD: 58

Base Curve: 4-6

## FEATURES:

Adjustable silicone nose pads, 0.50mm side shields included, versatile fit for all face types, lightweight metal frame

## COLORS:

Silver (SIL)

