



PROTECH
Medical

**PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

661/662S METALITE



MODEL: [RE-661/662S]

DESCRIPTION:

The model 557 is a metal frame featuring custom fitted 0.50mm lead strip side shields that create a seal to the front lenses for complete lateral coverage. Its gold color and retro style offer a classic look while non-slip, silicone nose pads provide easy adjustability and allow the frame to fit a wide variety of face types.

DESIGN:

Frame style: Metals

Front lens: 0.75mm Pb Glass

Frame Material: Optical Steel

Weight: 54g | 70g (662S)

Frame size: 52-18-144 | 52-18-145 (662S)

Side Shields: 0.50mm (included with 662S)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive

Range: +8.00 to -9.00

Minimum Progressive PD: 62

Base Curve: 4-6

FEATURES:

Adjustable silicone nose pads, 0.50mm side shields included, versatile fit for all face types, lightweight metal frame

COLORS:

Gold (GLD)

