



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

70 ASTROFLEX

MODEL: [RE-70S]



DESCRIPTION:

Impact resistance and flexibility come together in this high quality goggle. The frame features molded side shielding and is fitted to the brow for increased splash protection. The temple length and frame angle are adjustable for a uniform custom fit. Model 70S comes with leaded side shields for necessary lateral x-ray coverage.

DESIGN:

Frame style: Goggles
Front lens: 0.75mm Pb Glass
Frame Material: Plastic

Weight: 88g
Frame size: 54-19-Adjustable
Side Shields: 0.50mm Pb strips (optional)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -10.00

Minimum Progressive PD: 62
Base Curve: 4

FEATURES:

Wide field of view, Side Shields, Adjustable temples

