



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# 90 FITOVER

MODEL: [ RE-90 ]



## DESCRIPTION:

Designed to fit over prescription eyewear, model 90 fitovers are an excellent choice for those requiring protection against x-ray and liquid splashes. The large, clear side shields protect the user from lateral exposure while eliminating the "tunnel" effects of look-alikes.

## DESIGN:

Frame style: Fitover

Front lens: 0.75mm Pb Lead Glass

Frame Material: Plastic

Weight: 112g

Frame size: 66-14-120

Side Shields: 0.50mm (standard) or 0.75mm Pb Lead Glass

## PRESCRIPTION:

Options: N/A

Range: N/A

Minimum Progressive PD: N/A

Base Curve: 4

## FEATURES:

Large clear side shields

## COLORS:

Matte Black (MBLK)

