



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# 98 SUPERLITE

MODEL: [ RE-98 ]

## DESCRIPTION:

The 98 Superlite is an ultra light-weight frame designed for comfort. Its high curvature provides enhanced lateral/side protection against x-ray. Rubberized temples provide a secure, comfortable fit.

## DESIGN:

Frame style: Wraparound  
Front lens: 0.75mm Pb Glass  
Frame Material: Nylon

Weight: 63g  
Frame size: 63-15-120  
Side Shields: None

## PRESCRIPTION:

Options: None  
Range: Not applicable

Minimum Progressive PD: Not applicable  
Base Curve: 6-8

## FEATURES:

Rubberized nose pads and temples, lightweight nylon frame, highly curved lenses for optimal x-ray protection

## COLORS:

Brown (BRN)



BRN