



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

99 ULTRALITE



MODEL: [RE-99/99SM/99R]

DESCRIPTION:

The 99 Ultralites provide a wide area of protection against x-ray. Their high curvature provide adequate lateral protection which eliminates the need for side shields. The 99 comes in 2 sizes and an optional soft rubber nose bridge for additional support and comfort (99R). (Color blue not available in small size 99SM).

DESIGN:

Model: 99 / 99R (nose bridge)
Frame style: Wraparound
Front lens: 0.75mm Pb Glass
Frame Material: Plastic
Weight: 67g
Frame size: 64-18-140
Side Shields: N/A

Model: 99SM (small)
Frame style: Wraparound
Front lens: 0.75mm Pb Glass
Frame Material: Plastic
Weight: 59.5g
Frame size: 56-18-140
Side Shields: 0.50mm Pb strips (optional)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -10.00
Minimum Progressive PD: 64
Base Curve: 6

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -10.00
Minimum Progressive PD: 58
Base Curve: 6

FEATURES:

Rubber Nose Bridge (99R), Available in regular and small sizes(99SM), Rubberized Temples, Wide Field of View

COLORS:

Silver (SIL), Brown (BRN), Black (BLK), Blue (BLU), Fuschia (FUSC), Light Blue (LBLU), Light Pink (LPNK), Light Purple (LPUR), Teal (TEAL)

