



PROTECH
Medical

PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

ABACO DOCKSIDE



MODEL: [RE-AB-DOCK]

DESCRIPTION:

The Abaco Dockside lead glasses have a bold square style in an impact-resistant, polycarbonate frame featuring premium metal hinges. These fashion-forward and protective frames are both durable and lightweight.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Polycarbonate

Weight: 72g
Frame size: 58-19-135
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 60
Base Curve: 4-6

FEATURES:

Lightweight frame, versatile fit for all face types (unisex)

COLORS:

Black (BLK), Crystal Clear (CLR)



CLR



BLK