



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

ABACO JESSE

MODEL: [RE-AB-JES]

DESCRIPTION:

The Abaco Jesse lead glasses are unisex, universally flattering frames that are made of lightweight polycarbonate. These protective eyewear feature characteristic rivets at the temples and extremely flexible spring hinges for a snug, comfortable fit.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Polycarbonate

Weight: 65g
Frame size: 56-18-135
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 60
Base Curve: 4-6

FEATURES:

Lightweight frame, versatile fit for all face types (unisex)

COLORS:

Matte Black (MBLK)

