



**PROTECH**  
Medical

PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# ABACO KAI

MODEL: [ RE-AB-KAI ]



## DESCRIPTION:

Made of impact-resistant, lightweight polycarbonate, the Abaco Kai Lead glasses feature premium metal hinges and are molded around the ears to provide a comfortable fit. These unisex, protective frames complement and fit virtually all faces from small to large.

## DESIGN:

**Frame style:** Brand name - Classic fit  
**Front lens:** 0.75mm Pb Glass  
**Frame Material:** Polycarbonate

**Weight:** 65g  
**Frame size:** 56-18-135  
**Side Shields:** 0.50mm Strips (optional)

## PRESCRIPTION:

**Options:** Single Vision, Bifocal, Progressive  
**Range:** +8.00 to -9.00

**Minimum Progressive PD:** 60  
**Base Curve:** 4-6

## FEATURES:

Lightweight frame, versatile fit for all face types (unisex)

## COLORS:

Matte Black (MBLK), Crystal Clear (CLR)



CLR



MBLK