



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

ABACO TIKI

MODEL: [RE-AB-TIKI]



DESCRIPTION:

The Abaco Tiki is a classic-style frame made of a wood-textured, durable and lightweight polycarbonate. These lead glasses feature characteristic rivets at the temples, standard hinges, and are molded around the ears for added comfort. The Tiki protective eyewear are unisex and universally flattering to fit all faces from small to large.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Polycarbonate

Weight: 70g
Frame size: 54-21-140
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 60
Base Curve: 4-6

FEATURES:

Lightweight frame, versatile fit for all face types (unisex)

COLORS:

Dark wood (**WOOD**)

