



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

ABACO WAIKIKI

MODEL: [RE-AB-WAI]

DESCRIPTION:

The classic Abaco Waikiki lead glasses feature a lightweight polycarbonate frame to provide durability and a comfortable fit. The Abaco leaded eyewear feature characteristic rivet details on the front, premium metal hinges, and they are molded around the ears for added comfort. Waikiki glasses look good on virtually everyone.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Polycarbonate

Weight: 68g
Frame size: 54-21-140
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 60
Base Curve: 4-6

FEATURES:

Lightweight frame, versatile fit for all face types (unisex)

COLORS:

Tortoise (**TORT**), Shiny Black (**SBLK**)



TORT



SBLK