



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

AIRBORNE

MODEL: [RE-AIRB]

DESCRIPTION:

The Airborne is fully loaded with an Asian fit nose bridge, removable facial cavity, and curved lenses for optimal lateral x-ray protection.

DESIGN:

Frame style: Wraparounds

Front lens: 0.75mm Pb Glass

Frame Material: Polycarbonate

Weight: 74.7g

Frame size: 57-21-145

Side Shields: 0.50mm Pb Strips (optional)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive

Range: +8.00 to -10.00

Minimum Progressive PD: 62

Base Curve: 6

FEATURES:

High curve lenses for optimal x-ray protection, removable facial cavity foam insert, rubberized asian-fit nose bridge

COLORS:

Clear/Black (CLR)

