



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# CENTRIST

MODEL: [ RE-CENT ]

## DESCRIPTION:

The Centrist is a rectangular style safety frame made from high quality nylon and has permanent extended side shields for added x-ray protection. The sleek black with clear insert frame has an adjustable nosepiece that provides comfort and a secure fit.

## DESIGN:

Frame style: Wraparounds  
Front lens: 0.75mm Pb Glass  
Frame Material: Nylon

Weight: 88g  
Frame size: 55-17-140  
Side Shields: 0.50mm (standard) and optional extended side shields

## PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive  
Range: +8.00 to -10.00

Minimum Progressive PD: 62  
Base Curve: 6

## FEATURES:

Extended side shields, rubberized and adjustable asian-fit nose bridge

## COLORS:

Clear/Black (BLK)

