



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# CIRCUIT

MODEL: [ RE-CIRC ]

## DESCRIPTION:

The Circuit frame features wrap-around lenses which maximize lateral protection without the use of side shields. Rubberized temples and nose pads provide excellent comfort and grip. The Circuits modern style combined with its optimized functionality ensure long lasting quality and protection.

## DESIGN:

Frame style: Wraparound  
Front lens: 0.75mm Pb Glass  
Frame Material: Polycarbonate

Weight: 79g  
Frame size: 60-18-135  
Side Shields: 0.50mm Pb Strips (optional)

## PRESCRIPTION:

Options: None  
Range: Not applicable

Minimum Progressive PD: Not applicable  
Base Curve: 6

## FEATURES:

Rubberized nose pads and temples, versatile fit for all face types, wide field of view

## COLORS:

Silver/Black (SIL), Red/Gray (RED), Black/Gray (BLK)



SIL



BLK



RED