



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

CLASSIC METAL

MODEL: [RE-008]

DESCRIPTION:

Available in a beautiful light-blue color, this retro style frame is an aesthetically pleasing option for metal wire frames. The classic metal frame come standard with 0.50mm LE side shields for maximum lateral coverage. Non-slip, silicone nose pads offer easy adjustability and allow the frame to fit a wide variety of face types.

DESIGN:

Frame style: Metals
Front lens: 0.75mm Pb Glass
Frame Material: Optical Steel

Weight: 64g
Frame size: 54-19-140
Side Shields: 0.50mm (included)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 58
Base Curve: 4-6

FEATURES:

Adjustable silicone nose pads, 0.50mm side shields included, versatile fit for all face types, lightweight metal frame

COLORS:

Gold (GLD), Nano Blue (NBLU)

