



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

COMET

MODEL: [RE-COMET]



DESCRIPTION:

A sleek and versatile relative to the Proguard LeadR's frame, the Comet is an affordable option providing excellent coverage against x-ray. The Comets unique side shield configuration allows for lead shielding to be placed on the frame exterior and wrap around the entire lens area, providing a complete seal from scatter radiation.

DESIGN:

Frame style: Goggles
Front lens: 0.75mm Pb Glass
Frame Material: Plastic
Weight: 78.3g
Frame size: 59-19-114
Side Shields: 0.50mm Pb strips, Extended Side Shields (optional)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -10.00
Minimum Progressive PD: 64
Base Curve: 4

FEATURES:

Wide field of view, rubberized Asian fit nose bridge, rubberized temples, foam top bar, extended side shields

COLORS:

Blue/Black (BLU), Pink (PNK), Orange/Gray (ORNG)



BLU



PNK



ORNG