



PROTECH
Medical

**PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

FITGUARD ALPHA

MODEL: [RE-FGUARDA]

DESCRIPTION:

The Fitguard fitover is a versatile, high quality nylon frame. Featuring Flexfit technology, the fitover's temples can be benched (adjusted) easily for a secure, comfortable fit. The XL model is ideal for fitting over wide prescription eyewear while still providing a secure and comfortable fit. With the most balanced weight distribution, the XL is our most popular fitover model. Blacked out side shields provide a streamlined look as well as 0.50mm lead equivalent protection.

DESIGN:

Frame style: Fitover

Front lens: 0.75mm Pb Lead Glass

Frame Material: Plastic

Weight: 91g

Frame size: 69-17-145

Side Shields: 0.50mm (standard) Pb Strips

PRESCRIPTION:

Options: N/A

Range: N/A

Minimum Progressive PD: N/A

Base Curve: 4

FEATURES:

Large, blacked-out side shields

COLORS:

Matte Black (MBLK)

