



PROTECH
Medical

**PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

FITGUARD BETA



MODEL: [RE-FGUARDB]

DESCRIPTION:

The Fitguard fitover is a versatile, high quality nylon frame. Featuring Flexfit technology, the fitover's temples can be benched (adjusted) easily for a secure, comfortable fit. The Beta size is ideal for fitting over standard sized prescription eyewear with square or rounded lenses. Blacked out side shields provide a streamlined look as well as 0.50mm lead equivalent protection.

DESIGN:

Frame style: Fitover
Front lens: 0.75mm Pb Lead Glass
Frame Material: Plastic

Weight: 108g
Frame size: 64-18-150
Side Shields: 0.50mm (standard) Pb Strips

PRESCRIPTION:

Options: N/A
Range: N/A

Minimum Progressive PD: N/A
Base Curve: 4

FEATURES:

Large, blacked-out side shields

COLORS:

Matte Black (**MBLK**)

