



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

FULL FACE SHIELD



MODEL: [450FPM]

DESCRIPTION:

The Full Face Shield uses a lead impregnated protective acrylic shield suspended from a dual adjustable headpiece. The Full face shield provides maximum x-ray protection due to its larger lead impregnated shield. Because of the multi-adjustable headpiece, you can angle the mask to protect from radiation in many different positions as well as adjust the size to comfortably fit your head. *Mask should be worn with leaded eyewear protection.

DESIGN:

Model: 450FPM
Frame style: Face Shield
Front lens: 0.10mm Lead Acrylic

Frame Material: Plastic
Weight: 450g
Frame size: Adjustable Head Piece

PRESCRIPTION:

Options: N/A
Range: N/A

Minimum Progressive PD: N/A
Base Curve: N/A

