



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# MAKO

MODEL: [ RE-MAKO ]

## DESCRIPTION:

The Mako frame has highly curved lenses which maximize radiation protection coverage. Integrated side shields equipped with lead strip side shields provide essential lateral x-ray protection.

## DESIGN:

Frame style: Wraparound  
Front lens: 0.75mm Pb Glass  
Frame Material: Polycarbonate

Weight: 69.3g  
Frame size: 65-15-130  
Side Shields: 0.50mm Pb Strips

## PRESCRIPTION:

Options: None  
Range: Not applicable

Minimum Progressive PD: Not applicable  
Base Curve: 8

## FEATURES:

Rubberized nose pads and temples, highly curved frame for optimal x-ray protection, integrated lateral side shields

## COLORS:

Gray/Green (GRY)



GRY