



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

MALIBU

MODEL: [RE-MALIB]

DESCRIPTION:

The Malibu frame modernizes the classic wayfarer design and is offered in several tasteful color selections.

DESIGN:

Frame style: Classic

Front lens: 0.75mm Pb Glass

Frame Material: Polycarbonate

Weight: 68.1g

Frame size: 54-20-155

Side Shields: None

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive

Range: +8.00 to -10.00

Minimum Progressive PD: 58

Base Curve: 4

FEATURES:

Affordable and fashion forward

COLORS:

Matte Blue (MBLU), Matte Black (MBLK), Matte Brown (MBRN), Shiny Black (SBLK)



SBLK



MBLU



MBLK



MBRN