



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

MODERN METAL

MODEL: [RE-002]

DESCRIPTION:

The model 002 features a thin wire frame with contoured side shields. The ultra light-weight model 002 comes standard with 0.50mm LE side shields for maximum lateral coverage. Non-slip, silicone nose pads offer easy adjustability and allow the frame to fit a wide variety of face types.

DESIGN:

Frame style: Metals
Front lens: 0.75mm Pb Glass
Frame Material: Optical Steel

Weight: 73g
Frame size: 48-20-135
Side Shields: 0.50mm (included)

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 58
Base Curve: 4-6

FEATURES:

Adjustable silicone nose pads, 0.50mm side shields included, versatile fit for all face types

COLORS:

Gold (GLD), Light Gunmetal (GUNM)



GLD



GUNM

