



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone:	561.627.9769
Email(USA):	sales@protechmed.com
Email(Intl.):	global@protechmed.com
Email(LatAm.):	brazil@protechmed.com

# NIKE BRAZEN

MODEL: [ RE-N-BR ]

## DESCRIPTION:

The Nike Brazen lead glasses feature a ventilated nose bridge for grip, comfort, and reduced fogging. The Brazen's frame is constructed of lightweight nylon and offers superior coverage and protection.

## DESIGN:

**Frame style:** Brand name - Wraparound fit  
**Front lens:** 0.75mm Pb Glass  
**Frame Material:** Nylon

**Weight:** 75g  
**Frame size:** 60-16-140  
**Side Shields:** 0.50mm Pb strips (optional)

## PRESCRIPTION:

**Options:** Single Vision, Bifocal, Progressive  
**Range:** +8.00 to -9.00

**Minimum Progressive PD:** 64  
**Base Curve:** 6

## FEATURES:

Rubberized ventilated nose bridge, lightweight nylon frame, superior wraparound coverage for optimal protection

## COLORS:

Matte Black with blue (MBLK), Shiny Black (SBLK)



MBLK



SBLK