



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

NIKE MAVERICK

MODEL: [RE-N-MAV]

DESCRIPTION:

Utilizing Nike's unique nose bridge design, the Maverick lead glasses have reduced fogging and maximized comfort during long fluoroscopic procedures.

DESIGN:

Frame style: Brand name - Wraparound fit
Front lens: 0.75mm Pb Glass
Frame Material: Nylon

Weight: 70.3g
Frame size: 59-15-145
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 64
Base Curve: 6

FEATURES:

Unique nose bridge design for fog reduction, lightweight frame

COLORS:

Matte Black (MBLK)

