



PROTECH
Medical

**PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

NIKE PREMIER

MODEL: [RE-N-PRMR]

DESCRIPTION:

The Nike Premier lead glasses offer a diverse selection of colors and tasteful design elements. The Premier protective eyewear feature large rectangular lenses to provide adequate x-ray protection and coverage.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Nylon

Weight: 67.8g
Frame size: 60-13-135
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 64
Base Curve: 6

FEATURES:

Large rectangular lenses, lightweight frame

COLORS:

Matte Anthracite (**MANTH**), Matte Tortoise (**MTORT**), Matte Black (**MBLK**)



MANTH



MTORT



MBLK