



PROTECH
Medical

**PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL**

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

NIKE RABID

MODEL: [RE-N-RABID]

DESCRIPTION:

The Nike Rabid lead glasses feature a high-curvature sports frame providing superior wraparound coverage and protection. These protective eyewear boast a lightweight nylon and wide lenses for clarity at all angles during any procedure. Nike's Secure Wrap Technology features grip zones throughout the frame for a better overall fit, comfort, and stability.

DESIGN:

Frame style: Brand name - Wraparound fit
Front lens: 0.75mm Pb Glass
Frame Material: Nylon

Weight: 63g
Frame size: 63-14-145
Side Shields: Not available

PRESCRIPTION:

Options: None
Range: Not applicable

Minimum Progressive PD: Not applicable
Base Curve: 6

FEATURES:

High curvature for wraparound coverage, lightweight frame, superior grip for stability and comfort

COLORS:

Matte black with white swoosh (MBLK)

