



PROTECH LEADED EYEWEAR INC,  
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2  
Lake Park, FL 33403

Telephone: 561.627.9769  
Email(USA): sales@protechmed.com  
Email(Intl.): global@protechmed.com  
Email(LatAm.): brazil@protechmed.com

# OAKLEY CROSSLINK ZERO



MODEL: [ RE-OA-CROSS ]

## DESCRIPTION:

The Crosslink Zero has fixed temples with a three-Point Fit to ensure lenses are aligned for visual correction. Wire cores hidden inside the temples of this lightweight frame offer adjustability for a perfect fit.

## DESIGN:

**Frame style:** Brand name - Wraparound fit  
**Front lens:** 0.75mm Pb Glass  
**Frame Material:** Durable O Matter™

**Weight:** 61.6g  
**Frame size:** 58-16-138  
**Side Shields:** Not available

## PRESCRIPTION:

**Options:** Single Vision, Bifocal, Progressive  
**Range:** +8.00 to -9.00

**Minimum Progressive PD:** 56  
**Base Curve:** 4-6

## FEATURES:

Rubberized ventilated nose bridge, lightweight nylon frame, superior wraparound coverage for optimal protection

## COLORS:

Polished Gray Shadow (**GRY**), Matte Black (**MBLK**)



GRY



MBLK