



PROTECH
Medical

PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

OAKLEY FIVE SQUARED

MODEL: [RE-OA-FIVESQ]

DESCRIPTION:

The Oakley Five Squared leaded glasses use Condensed Cranial Geometry for a tapered architecture that feels as good as it looks. Dimensional reliefs, metal icons and sculpturally integrated hinges take the whole idea of street style and amp it with attitude.

DESIGN:

Frame style: Brand name - Wraparound fit
Front lens: 0.75mm Pb Glass
Frame Material: Durable O Matter™

Weight: 71g
Frame size: 54-20-135
Side Shields: Not available

PRESCRIPTION:

Options: None
Range: Not applicable

Minimum Progressive PD: Not applicable
Base Curve: 8

FEATURES:

Wraparound coverage for optimal protection

COLORS:

Gray Smoke (GRY), Polished Black (PBLK)



PBLK

OAKLEY



GRY