



PROTECH
Medical

PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

OAKLEY HOLBROOK

MODEL: [RE-OA-HOL]

DESCRIPTION:

The Oakley Holbrook lead glasses combine traditional design with a modern flare.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Durable O Matter™

Weight: 66.5g
Frame size: 56-18-135
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 60
Base Curve: 4-6

FEATURES:

Lightweight frame

COLORS:

Polished Clear (CLR), Matte Brown Tortoise (MTORT), Satin Black (SBLK)



CLR



SBLK



MTORT