



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

OAKLEY LITEBEAM

MODEL: [RE-OA-LITE]

DESCRIPTION:

Smooth details and streamlined design provide a timeless look with the Oakley Litebeam lead glasses. Composed of lightweight O Matter™ material, the Litebeam protective eyewear feature adjustable temples and come with several nosebridge sizes to ensure maximum comfort.

DESIGN:

Frame style: Brand name - Classic fit
Front lens: 0.75mm Pb Glass
Frame Material: Durable O Matter™

Weight: 59.6g
Frame size: 55-18-140
Side Shields: Not available

PRESCRIPTION:

Options: Single Vision, Bifocal, Progressive
Range: +8.00 to -9.00

Minimum Progressive PD: 60
Base Curve: 4

FEATURES:

Adjustable temples for perfect fit, Lightweight frame

COLORS:

Satin Black (SBLK)

