



PROTECH LEADED EYEWEAR INC,
DBA PROTECH MEDICAL

1360 North Killian Drive, Unit 2
Lake Park, FL 33403

Telephone: 561.627.9769
Email(USA): sales@protechmed.com
Email(Intl.): global@protechmed.com
Email(LatAm.): brazil@protechmed.com

OAKLEY STRAIGHTLINK

MODEL: [RE-OA-STRGHT]

DESCRIPTION:

The Oakley Straightlink lead glasses feature rubberized earsocks and nosepads for a no-slip grip with perspiration, making this frame great for long procedures. These protective glasses have a Three-Point Fit which holds lenses in precise optical alignment while eliminating pressure points.

DESIGN:

Frame style: Brand name - Wraparound fit
Front lens: 0.75mm Pb Glass
Frame Material: Durable O Matter™

Weight: 77.2g
Frame size: 61-19-130
Side Shields: Not available

PRESCRIPTION:

Options: None
Range: Not applicable

Minimum Progressive PD: Not applicable
Base Curve: 8

FEATURES:

Rubberized temples and nose pads for comfort and grip, optimized weight distribution and alignment, lightweight frame

COLORS:

Matte Black (MBLK)

